

ASMC 2016



Measurement of Nanoparticles on Parts: Figures of Merit of a Liquid-borne Particle Counter

Bin Liu and Shi Liu

ChemTrace Analytical Services, Fremont, California, USA

Abstract

- Measuring surface particles extracted from parts by liquid-borne particle counter (LPC) effective in evaluating surface cleanliness of semiconductor tool components (parts)
- Extended previous framework for evaluating commercial LPCs to a LPC counting particles smaller than 200 nm with close to 65% counting efficiency
- Evaluated counting accuracy, precision, stability, and limit of detection, among others, of a commercial LPC
- Counting efficiency not considered, counting accuracy <70% from results based on PSL count standard
- After counting efficiency correction, counting error <3%. Measured imprecision ≤3%. The linear dynamic range from 15 to 1,000 particles/mL.

Results




Fig. 1. Counting accuracy evaluated by measuring a PSL count standard. Ten measurements of count standard (0.231µm nominal size, certified value 1,002 particles/mL) average 642 particles/mL. Without counting efficiency correction, accuracy based on this result is 64%; after correction, accuracy is 98.6%.




Fig. 2. Linear dynamic range determined with a set of dilutions made from a quartz extract. The calculated concentrations range from 400 to 1,200 particles/mL. Measured linear dynamic range from 15 to approximately 1,000 particles/mL. Measured concentrations above this range start to deviate.




Fig. 3. Measurements showing high precision of this counter. 25 repeated measurements in Experiment 1 and 35 repeated measurements in Experiment 2. Both experiments show ≤3% imprecision.




Fig. 4. Extracts of different materials measured for stability. Water and yttria extracts were relatively stable; SiC extract stable after Day 2; Alumina extract different everyday.

Table I: Spike recovery experiments performed by spiking a water sample with different extracts, respectively. Both spiked and unspiked samples were measured. Spike recoveries range from 92% to 119% for particle concentrations above 15 particles/mL (LOQ).

Material	Spiking Level (particles/mL) / Recovery (%)				
	≥ 0.1 µm	≥ 0.15 µm	≥ 0.2 µm	≥ 0.3 µm	≥ 0.5 µm
Quartz	884	298	129	25	6
	99%	99%	98%	92%	67%
Yttria	1064	475	235	51	12
	102%	103%	105%	107%	101%
Al ₂ O ₃ Ceramic	1107	443	208	44	5
	101%	105%	108%	119%	139%
Doped Yttria	757	452	293	121	41
	108%	107%	106%	109%	112%
Aluminum	954	494	286	99	30
	108%	110%	110%	108%	98%




Fig. 5. Gage R&R performed per EMP III with ten different extracts. Extracts measured by two operators in triplicate for a total of 60 measurements. According to EMP III, the measurement process is a “first-class monitor”

Conclusion

- Counting efficiency affects measurement accuracy
- Linear dynamic range different from manufacturer spec of maximum 1,200 particles/mL
- Evaluating surface particles extracted from parts by LPC is effective only if LPC is fit for particle contamination measurement.
- Each LPC might have a different counting efficiency, direct comparison of uncompensated result is meaningless.
- The effect of counting efficiency on different counters and different materials will be studied in further details in future work.

References

T. Hattori (Ed.), Ultraclean Surface Processing of Silicon Wafers (Springer-Verlag, Berlin Heidelberg, 1998)
 K.L. Mittal, R. Jaiswal (Ed.), Particle Adhesion and Removal (Scrivener Publishing LLC, Massachusetts, 2015)
 S. H. Tan, N. Chen, S. Liu, K. Wang, Contamination-Free Manufacturing: Tool Component Qualification, Verification and Correlation with Wafers (Characterization and Metrology for ULSI Technology, 2003)
 R. P. Donovan (Ed.), Particle Control for Semiconductor Manufacturing (Marcel Dekker, Inc., New York, 1990).
 K. Kondo, Latest Technology and Standardization Trends for Liquid-borne Particle Counter (pp 449, Bunseki, Vol. 9, 2012)
 K. Kondo, Measurement and Control of Liquid-borne Particles in the Semiconductor Manufacturing Process (p. 37, Cleaning Technology, Vol. 22, No. 12)
 S. Liu, B. Liu, Particle Measurement with a Liquid-Borne Particle Counter: Analytical Figures of Merit (pp 157 – 160, Solid State Phenomena, Vol. 219, 2015)
 Automotive Industry Action Group, Measurement Systems Analysis; Reference Manual (4th ed.). (Automotive Industry Action Group, Detroit, 2010)
 Donald J. Wheeler, EMP III: Evaluating the Measurement Process & Using Imperfect Data (SPC Press, Tennessee, 2006)
 ISO 21501-2:2007(E). Determination of particle size distribution – Single particle light interaction methods – Part 2: Light scattering liquid-borne particle counter (2007).